# optimis Fusion

NEXT GENERATION SLT/YAG LASER



**Integrated SLT/YAG Platform** 





Quantel Medical's cutting edge laser technology has brought constant and multiple innovations to ophthalmologists worldwide **since 1993**.

Optimis Fusion combines traditional YAG photodisruption treatments and advanced SLT photo-regenaration therapy in one integrated platform.

# Optimis Fusion from Quantel The Next Generation of SLT/YAG Technology

Quantel's proprietary technology is based on an engineering breakthrough, which enables infrared and visible green nanosecond pulses to be produced by a single dual-mode cavity reducing complexity and resulting in a more elegant, advanced combination laser platform.

## **Premium Optics**

- Optimized for viewing the anterior and posterior segment
- Convergent optical path
- Blue-free illumination for improved visibility

### New, Advanced Technology

- Highly efficient, powerful and reliable
- Dual-mode laser cavity:
   Optimized to produce two wavelengths from one cavity

## New, Advanced Design

- Integrated Platform
- Easy switch from SLT to YAG mode
- Small footprint for two powerful technologies



# **SLT Mode**



# Gentle Laser Delivery

#### No Thermal Damage. Repeatable.

Selective Laser Trabeculoplasty (SLT) targets melanin rich cells which absorb the laser energy to induce a healing response that begins to restore the function of the trabecular meshwork.

The laser energy is delivered in short, fixed pulses over a large, homogeneous spot size, resulting in no thermal damage. Due to the gentle delivery and no permanent destruction to the trabecular meshwork, SLT leaves the structures intact and can be repeated if necessary.

# **YAG Mode**

The YAG mode delivers high-performance photodisruption for capsulotomy and peripheral iridotomy. The gaussian laser beam profile allows for precise delivery, only requiring minimum energy levels, which avoids side effects such as lens pitting.

Precise tissue targeting is achieved with the advanced two-point aiming beam, allowing for fast and accurate targeting of the capsule; and the adjustable anterior and posterior offset feature, allowing for easy focus adjustment in front, behind or at the point of disruption.

# Precise and Powerful







# SLT

First choice treatment for managing glaucoma and lowering IOP.

Repeatable Proven Gentle and non-invasive Primary or adjunct therapy No compliance issues No side effects\* \*Compared to drug therapy

# optimis Fusion





# TECHNICAL SPECIFICATIONS

OPTIMIS FUSION

	<u> </u>	
Laser source:	Q-switched, frequency doubled Nd:YAG	Q-switched Nd:YAG
Wavelength:	532 nm	1064 nm
Energy levels/Power:	0.3 to 2 mJ	0.3 to 10 mJ (1 pulse)
Pulse duration/Exposure:	4 ns	4 ns
Spot Size:	400 microns	10 microns (86.4% energy in diameter 8 microns FWHM
Repetition rate:	2.5 Hz	2 Hz
Offset: (Anterior and Posterior)	N/A	-150μm, +30μm, +150μm
Pulse mode:	Single pulse	1, 2 and 3 pulses per shot

#### **COMMON TO BOTH**

Aiming beam: Red 635 nm, adjustable intensity

Slit lamp: CSO 980 - LED

Magnification: 3 positions: 10x, 16x, 25x

Cooling: Air cooled

Size (laser only): 64.7 (H) x 54 (W) x 48.5 (D) cm

25.5" (H) x 21.3" (W) x 19" (D)

Weight (laser only): 24.15 kg 53.1 lbs

53.1 lbs

Power requirements: 100-240 VAC, 250 VA, 50/60 Hz

Delivery Systems: N/A

#### **AVAILABLE MODELS**

Optimis Fusion system

**Poland** 

Quantel Medical

ul. Racławicka 93

Tel: +48 22 5210111

02-634 Warszawa - POLAND

Email: info@quantel-medical.pl

Optimis Fusion and VITRA 2 SingleSpot laser platform

Optimis Fusion and VITRA 2 MultiSpot laser platform

Specifications are subject to change without notice.

© 2020. Quantel Medical, Optimis Fusion is a registered trademark of Quantel Medical. All rights reserved.

## www.quantel-medical.com





#### Headquarters

Quantel Medical 11, rue du Bois Joli - CS40015 63808 Cournon d'Auvergne - FRANCE Tel: +33 (0)4 73 745 745 Email: contact@quantel-medical.fr

#### North America

Quantel Medical USA 2313 Ridge Road Suite 105A Rockwall, TX 75087 – USA Tel: +1 888 660 6726 Email: info@quantelmedical.com





BY LUMIBIRD

ISO 9001 : 2015 - ISO 13485 : 2016

**Representative Offices** Thailand, Chiang Mai Brazil, Rio De Janeiro